

# **Automate measurements with Python**

Martin Mihaylov

28. September 2023

# Inhaltsverzeichnis

<b>1</b>	<b>Introduction to the Automatic Measurements Python Libs</b>	<b>3</b>
1.1	Library and Code Structures . . . . .	3
1.2	Set up Python . . . . .	6
<b>2</b>	<b>How to connect an instrument to PC/Laptop</b>	<b>8</b>
<b>3</b>	<b>Script Structure</b>	<b>14</b>
<b>4</b>	<b>Example</b>	<b>28</b>

# 1 Introduction to the Automatic Measurements Python Libs

Automatic Measurements is a Python library dedicated to remotely controlling laboratory measurement instruments in two separate labs. The library is an ongoing project, continuously updated by me to address errors, malfunctions, and spelling mistakes as they are discovered. The library is hosted on the Uni Paderborn Git Server. To enhance functionality while ensuring code cleanliness and minimizing errors, I am currently the sole contributor responsible for submitting changes to the Git repository. However, I encourage everyone to submit suggestions and code improvements to streamline and lighten the codebase. You might wonder why Python is chosen over other programming languages like LabVIEW or MATLAB. Python leverages NiVisa, Serial Port, FPGA Control, and Socket libraries, which are also used by other programs. The advantage of Python lies in its freedom from license checkouts and package concerns. Most libraries used are free, except those provided by specific manufacturers like CST or Lumerical for EM simulations. While Python may require more lines of code, it offers a robust and license-free solution.

## 1.1 Library and Code Structures

The Automatic Measurements Library is organized into multiple folders and scripts, designed to facilitate correction, extension, and maintenance. Within each section of the code, primary functions are accompanied by informative docstrings, providing clear guidance to users. As an example, in Figure 1.1, you can observe a sample docstring.

The `main.py` file serves as the central hub for interconnecting all the libraries. Users can access and control all the instruments and functions through this file, which can be located directly in the Auto\_Measurement folder within the Git repository. Here is a list of the instruments that can be controlled by the Python script:

1. Anritsu spectrum Analyzer MS2760A
2. Anritsu Signal Generator MG3694C
3. Anritsu Vectro Analyzer MS4647B
4. Power Meter ThorLabs PM100D
5. Novoptel Laser LU1000
6. Yokogawa Optical Spectrum Analyzer AQ6370D

```

def Laserbeating_RF_Response_CoBrite(Instrument, SA_f_min, SA_f_max, Laser_f_opt, Laser_f_max, Step, resBW, LaserChannel, LaserPower, path, name, SA_...
...


Parameters
-----
Instrument : List
    List of Instruments
SA_f_min : float
    Spectrum Analyzer minimal frequency in Hz
SA_f_max : float
    Spectrum Analyzer minimal frequency in Hz
Laser_f_opt : float
    Optical frequency of the laser in THz. For Example:
        Laser_f_opt = 192.3421
Laser_f_max : float
    Maximal frequency of the laser in THz. For Example:
        Laser_f_max = 192.3421
Step : int
    Steps from Laser_f_opt to Laser_f_max. For Example:
        Steps = 2 , Laser_f_opt = 192.3000 Laser_f_max = 192.3400 ->
        -(Laser_f_max- Laser_f_opt)/Steps
resBW : float
    Resolution bandwidth in Hz
LaserChannel : int
    Laser channel selected. It can be only 1 and 2 for CoBrite Laser Unit
LaserPower : int
    Laser output power in dB.
SA_TraceNum : int
    Spectrum Analyzer Active Trace. Per Default = 1
path : str
    Where to save the file.
name : str
    Name of the file. The name is Combination of the function name and
    name that you give.
    For Example: Coupling_Stability_Test.csv

Raises
-----
ValueError
    DESCRIPTION.

Returns
-----
None data will be returned. A TXT File with the Power Detector data will
be created, CSV file with data will be created and SVG plot
of the measurement will be created and saved in the given in 'path' folder.

```

*Abbildung 1.1: Example of an docstring*

7. KEITHLEY Source Meter 2612
8. Power Supply KA3005 (All Programmable KA Series Power supplies !)
9. CoBrite Tunable Laser
10. AnaPico AG, APPH20G
11. 4-Channels Power Supply GPP4323

In the InstrumentControl folder, you'll find scripts for controlling various instruments. These scripts detail which functions of the instruments are included and which are not. Additionally, within the Auto\_Measurement directory, you'll discover separate folders for each instrument, each containing a control script and a test main.py file. These main.py files were used by me when extending the functionality of the given instrument. Furthermore, you'll find PDF files containing documentation for each instrument's functions. To identify which functions are already programmed for each instrument, navigate to Auto\_Measurement/Instruments Docu/\_build/html/index.html. Within this directory, you can select the instrument of your choice to view the list of included functions. An example from the MS2760A can be seen in Figure 1.2.

As you can see the instrument functions are mostly separated in two types:

1. Ask With ask.... The user can ask the instrument for a value or data. For example: `ask_Volt(1)` Will return the Voltage measured by an GPP4323 4-Channel Power Supply

```

set_Continuous(status)
    Parameters: status (str/int) – Stop/start Sweep. Can be ['ON', 'OFF', 1, 0]
    Raises:      ValueError – Error message
    Returns:
    Return type: None.

set_ContinuousMeas(state)
    Parameters: state (str/int) – Title: Sweep Type Description: Specifies
                whether the sweep/measurement is triggered continuously. If
                the value is set to ON or 1, another sweep/measurement is
                triggered as soon as the current one completes. If continuous
                is set to OFF or 0, the instrument remains initiated until the
                current sweep/measurement completes, then enters the 'idle'
                state and waits for the :INITiate[:IMMediate] command or for
                :INITiate:CONTinuous ON.
    Raises:      ValueError – Error message
    Returns:
    Return type: None.

set_DataFormat(status)
    Parameters: status (str) – Set Data Format status =
                ['ASCII', 'INTeger', 'REAL']
    Raises:      ValueError – Error message
    Returns:
    Return type: None.

set_DataPointCount(value)
    Parameters: value (int) –
        Title: Display Point Count
        Description: Changes the number of display points the
                    instrument currently measures. Increasing the number of
                    display points can improve the resolution of
                    measurements but will also increase sweep time.
        Default Value: 501 Range: 10 to 10001
    Raises:      ValueError – Error message
    Returns:
    Return type: None.

```

*Abbildung 1.2: Small cut from the index.html*

2. Set With set.... The user can set an Value or some parameter on the choosen instrument. For Example: `set_VoltageLimit(A, 2.34)` Will set the Voltage on Channel A to 2.34 for Keithley 2612 Power Meter.

Indeed, in addition to the standard control functions for the instruments, there are specialized functions with naming conventions like get.... or save..... These functions are designed to extract specific data and either save it to a designated directory on your computer or make it available in your Python Variable Explorer. It's important to thoroughly review the index.html documentation to gain a comprehensive understanding of all available functions for each instrument. This will enable you to leverage the full capabilities of the Automatic Measurements Library and efficiently manage data extraction and storage.

If you have any specific questions or need further assistance with any of these functions

or instruments, please feel free to ask.

## 1.2 Set up Python

Python is readily installed on the laboratory PCs in the Optic Laboratory, complete with all the necessary libraries. To access Python, simply type Spyder into the Windows search menu. Once Spyder loads, open the main.py file and follow the instructions provided in the Python terminal. Alternatively, you can refer to Chapter 4 for guidance. Spyder is the chosen Python Integrated Development and Learning Environment (IDLE) in our laboratory. While other options like PyCharm, Vim, Visual Studio Code, and more are available, Spyder is favored for its lightweight and user-friendly nature. It can be installed as part of the Anaconda package or separately. Spyder is well-suited for scientific projects and is relatively easy to learn, particularly if you are familiar with Octave or Matlab IDLE. If Python is not installed on your PC or laptop, you can download it from Python. Ensure that you download a **STABLE** Python release. Following this, you can download Spyder or Anaconda from Spyder and Anaconda. You may also explore other Python APIs for your convenience.

After downloading and installing Python, you'll need to install several libraries to use the scripts. The installation process for libraries can vary depending on whether you are using Anaconda or if you installed Python first and then Spyder IDLE. If you're using Windows or Linux Terminal, you can install libraries using the following command **pip install [Library Name]**. Most libraries can be found on PyPi or by conducting a Google search. If you're using Anaconda, you can install libraries either through the Anaconda Terminal or directly from the built-in Anaconda Library manager. Here is a list of Python libraries that were installed and used on the laboratory PC with Python version 3.6:

1. tk=8.6.10=he774522\_0
2. pandas=1.1.3=py38ha925a31\_0
3. numpy=1.19.2=py38hadc3359\_0
4. matplotlib=3.3.2=0
5. scipy=1.5.2=py38h14eb087\_0
6. seaborn=0.11.0=py\_0
7. sockets=1.0.0=pypi\_0
8. python-vxi11=0.9=pypi\_0
9. pyusb=1.2.1=pypi\_0
10. libusb=1.0.24b1=pypi\_0
11. pyvisa=1.11.3=pypi\_0

12. pyvisa-py=0.5.2=pypi\_0
13. ftd2xx=1.1.2=pypi\_0
14. ftd3xx-py3k=1.0=pypi\_0
15. oct2py=5.3.0=pypi\_0
16. octave-kernel=0.33.1=pypi\_0

To utilize the scripts effectively, it is essential to install these libraries. While not all of them are required every time you use the scripts, having them all installed from the beginning can save you the hassle of installing new libraries each time you work with a different instrument. The numbers that follow the names of each library in the list are release and version numbers, which you can disregard. You should have at least these versions of the libraries or newer. Please be aware that this manual was written in 2023, so it's possible that some of these libraries may not be supported or could undergo changes 10 or 20 years from now. In such cases, consulting online resources and documentation or seeking assistance through online communities will be valuable for resolving any future confusion.

## 2 How to connect an instrument to PC/Laptop

In this chapter, we will provide detailed information about the laboratory instruments and the required cables for connecting them, complete with accompanying pictures for clarity.



(a) Anritsu spectrum Analyzer MS2760A picture from above



(b) Anritsu spectrum Analyzer MS2760A USB-C cable for PC connection

*Abbildung 2.1:* Anritsu spectrum Analyzer MS2760A picture from above 2.1a and behind 2.1b



(a) Anritsu Signal Generator MG3694C




(b) Anritsu Signal Generator MG3694C uses Ethernet cable to connect to PC

*Abbildung 2.2:* Anritsu Signal Generator MG3694C picture from above 2.2a and behind 2.2b



(a) Anrtisu Vectro Analyzer MS4647B



(b) Anrtisu Vectro Analyzer MS4647B uses Ethernet cable to connec to PC

**Abbildung 2.3:** Anrtisu spectrum Analyzer MS2760A picture from above 2.3a and behind 2.3b



(a) Power Meter ThorLabs PM100D



(b) Power Meter ThorLabs PM100D picture from the side, this instruments uses mini USB to connect to PC

**Abbildung 2.4:** Power Meter ThorLabs PM100D picture from above 2.4a and the side 2.4b




(a) Novoptel Laser LU1000



(b) Novoptel Laser LU1000 uses USB-A to USB-B Cable

**Abbildung 2.5:** Novoptel Laser LU1000 front 2.5a and behind 2.5b instrument picture




(a) Yokogawa Optical Spectrum Analyzer AQ6370D




(b) Yokogawa Optical Spectrum Analyzer AQ6370D uses Ethernet cable to connect to PC

**Abbildung 2.6:** Anritsu spectrum Analyzer MS2760A front 2.3a and behind 2.3b pictures




(a) KEITHLEY Source Meter 2612




(b) KEITHLEY Source Meter 2612 uses RS232 to USB cable to connect to PC

*Abbildung 2.7: Anrtisu spectrum Analyzer MS2760A front 2.7a and behind 2.7b pictures*



(a) Power Supply KA3005



(b) Power Supply KA3005 uses USB-A to USB-B cable to connect to PC

*Abbildung 2.8: Power Supply KA3005 front 2.7a and behind 2.7b pictures*



(a) CoBrite Tunable Laser



(b) CoBrite Tunable Laser uses RS232 to USB cable to connect to PC

**Abbildung 2.9:** CoBrite Tunable Laser front 2.9a and behind 2.9b pictures




(a) AnaPico AG APPH20G




(b) AnaPico AG APPH20G uses Ethernet cable to connect to PC

**Abbildung 2.10:** AnaPico AG APPH20G front 2.10a and behind 2.10b pictures



(a) 4-Channels Power Supply GPP4323



(b) 4-Channels Power Supply GPP4323 uses  
USB-A to USB-B cable to connect to PC

*Abbildung 2.11*

# 3 Script Structure

In this chapter the script structure will be presented. Certainly, providing a clear script structure is a valuable resource for users. It will help users navigate and understand the scripts more effectively, making it easier for them to find the information they need.

The scripts can be started in the with the main.py file. There all the needed librarys will be at first loaded. After starting the main.py the user will be askt in the python terminal to give the number of instruments that he/she want to use. For example if i have an measurments set and i need to use the power meter and voltage source i will put 2. After thus an list with all the instruments will be presented and the user need to choose witch are the two instruments. All of the instruments are numerated and the user can see the names and numbers in the python terminal. An example of the code is given in the snip below:

```
1 import tkinter as tk
2 from tkinter import filedialog
3 from InstrumentControl.InstrumentSelect import InstInit
4 from InstrumentControl.FunctionCall import MeasFunctions
5
6
7
8
9 #
=====
10 # Select Instruments and Load Instrument Libraries
11 #
=====

12
13 print('''
14     ##### How many instruments will be used #####
15
16     ''')
17 numInstr = int(input('Number Instruments = '))
18 print('''
19
20     ##### How many instruments will be used #####
21
22     ''')

23
24 listInstruments = []
25
26
27 print('''
28     ##### Select which instumnets will be used #####
```

```

29
30     You will be ask to give the name of the instruments
31     manualy. For example Instrument name: MS4647B for
32     Anrtisu Signal Analyzer.
33
34
35     Table of instruments
36     Anritsu spectrum Analyzer MS2760A = 1
37     Anritsu Signal Generator MG3694C = 2
38     Anritsu Vectro Analyzer MS4647B = 3
39     Power Meter ThorLabs PM100D = 4
40     Novoptel Laser LU1000 = 5
41     Yokogawa Optical Spectrum Analyzer AQ6370D = 6
42     KEITHLEY Source Meter 2612 = 7
43     Power Supply KA3005 = 8
44     CoBrite Tunable Laser = 9
45     AnaPico AG,APPH20G = 10
46     4-Channels Power Suppy GPP4323 = 11
47
48
49
50
51 ##### Select which instumnets will be used #####
52     ')
53 instrumenNameList = [ '1' , '2' , '3' , '4' , '5' , '6' , '7' , '8' , '9' , '10' , '11' ]
54 Instrument = []
55 for i in range(numInstr):
56     instrumenName = input( 'Instrument name: ' )
57     if instrumenName in instrumenNameList:
58         listInstruments.append(instrumenName)
59     else:
60         raise ValueError( 'Invalid Instrument Selected' )
61
62
63 for elements in listInstruments:
64     if elements == '1':
65         SA = InstInit(elements)
66         Instrument.append(SA)
67         print( 'Anrtisu Spectrum Analyzer MS2760A is connected as SA' )
68
69     elif elements == '2':
70         SG = InstInit(elements)
71         Instrument.append(SG)
72         print( 'Anrtisu Signal Generator MG3694C is connected as SG' )
73
74     elif elements == '3':
75         VNA = InstInit(elements)
76         Instrument.append(VNA)
77         print( 'Anrtisu Vectro Network Analyzer MS4647B is connected as VNA' )
78
79     elif elements == '4':
80         PM = InstInit(elements)
81         Instrument.append(PM)

```

```

82     print('Power Meter ThorLabs PM100D is connected as PM')
83
84 elif elements == '5':
85     LU = InstInit(elements)
86     Instrument.append(LU)
87     print('Novoptel Laser LU1000 is connected as LU')
88
89 elif elements == '6':
90     OSA = InstInit(elements)
91     Instrument.append(OSA)
92     print('Yokogawa AQ6370D is connected as OSA')
93
94 elif elements == '7':
95     KA = InstInit(elements)
96     Instrument.append(KA)
97     print('KEITHLEY Source Meter 2612 is connected as KA')
98
99 elif elements == '8':
100    SerialNum = [ 'KORAD KA3005P V5.8 SN:03379314' , 'KORAD KA3005P V5.8
101 SN:03379289' , 'RND 320-KA3005P V2.0' ]
102    PS = InstInit(elements)
103    data = PS.getIdn()
104    print(data)
105    if data in SerialNum:
106        if data == SerialNum[0]:
107            PS1 = PS
108            print('Power Supply KA3005P is connected as PS1')
109            Instrument.append(PS1)
110        elif data == SerialNum[1]:
111            PS2 = PS
112            print('Power Supply KA3005P is connected as PS2')
113            Instrument.append(PS2)
114        if data == SerialNum[2]:
115            PS3 = PS
116            print('Power Supply KA3005 is connected as PS3')
117            Instrument.append(PS3)
118
119 elif elements == '9':
120     CO = InstInit(elements)
121     Instrument.append(CO)
122     print('CoBrite Tunable Laser is connected as CO')
123 elif elements == '10':
124     AP = InstInit(elements)
125     Instrument.append(AP)
126     print('AnaPico AG,APPH20G is connected as AP')
127 elif elements == '11':
128     GPP = InstInit(elements)
129     Instrument.append(GPP)
130     print('GW Instek,GPP4323 is connected as GPP')
131 else:
132     raise ValueError('Invalid Instrument Selected')
133
134

```

```

135 #
136 # Select Path for Save Data (pop up window)
137 #
138 root = tk.Tk()
139 path = filedialog.askdirectory(parent = root)
140 root.destroy()

```

The constant prompting for commands in the terminal serves as a practical solution to the challenge of not having a graphical user interface. Once the user specifies the number of instruments they wish to use and selects the instrument(s) by name and number, the main.py script will invoke InstrumentSelect.py. Within InstrumentSelect.py, an attempt is made to establish connections with the selected instrument(s). Each instrument may require a different form of connection, such as Ethernet, USB, etc. Therefore, the script begins by inspecting the available PC ports to locate the desired instrument. To streamline the usability of the code, instruments utilizing USB connections are identified by their serial numbers. This is crucial because when a new instrument is bought for one of the laboratories, its serial number must be added to the system for user-friendly connectivity. During the process, informative messages will be displayed in the Python terminal to indicate that the system is attempting to establish connections. It's worth noting that these messages may appear multiple times if the PC has multiple USB ports. As an example, the following code snippet illustrates how power supplies may be connected within the InstrumentSelect.py script:

```

1 def PowerSupply():
2     from InstrumentControl.RD3005 import RD3005
3     from InstrumentControl.KA3005 import KA3005
4     from InstrumentControl.KA3005p import KA3005p
5
6     SerialNum = [ 'KORAD KA3005P V5.8 SN:03379314' , 'KORAD KA3005P V5.8 SN
7 :03379289' , 'RND 320-KA3005P V2.0' ]
8     #Prnt all instruments connected to the COM-Ports .
9     #Needed to set later
10    import serial.tools.list_ports
11    ports = serial.tools.list_ports.comports()
12    COM_List = []
13    Port_ = None
14    for port, desc, hwid in sorted(ports):
15        # print("{}: {} [{}]").format(port, desc, hwid))
16        COM_List.append(port)
17
18    PowerInstr = 0
19    for data in list(COM_List):
20        while PowerInstr == 0:
21            try:
22                PS = RD3005(data)
23                PowerInstr = PS.getIdn().split("\n")[0]
24                if PowerInstr in SerialNum:

```

```

25         PowerInstr = 1
26         Port_ = data
27         break
28     else:
29         PowerInstr = 0
30         print("Scanning COM Ports for Instrument !")
31     except serial.SerialException as e:
32 #There is no new data from serial port
33         print("Scanning COM Ports for Instrument !")
34     except TypeError as e:
35         #Disconnect of USB->UART occured
36         print("Scanning COM Ports for Instrument !")
37     except visa.VisaIOError as e:
38         print("Scanning COM Ports for Instrument !")
39     except AttributeError:
40         pass
41     else:
42         break
43     break
44 CheckInstrName = None
45 CheckInstrName = PS.getIdn().split("\n")[0]
46 PS.Close()

```

After successfully connecting the instrument, it is passed back to the main.py script. In the main.py terminal, you will receive a notification like: **Power Supply KA3005P is connected as PS1**. This notification confirms the successful connection of, for example, the Power Supply KA3005P, which is identified as 'PS1.' Once the connection process for all selected instruments is completed, a final dialog window will appear, prompting you to specify the directory where you intend to save measurement data. Subsequently, you can easily interact with the instruments by using commands such as: **PS1.Instrument function**. In addition, the Python scripts incorporate a secondary section where predefined measurement functions are available. These functions can be initiated with:

- 1 MeasFunctions(Instrument ,path)

In this section, we provide a detailed list of all the pre-coded functions available for your convenience. After selecting a function in the Python terminal, you will be prompted to provide the initial parameters required to perform measurements. You can refer to Chapter 4 for a practical demonstration of how this process works. The functions are selected in FunctionCall.py, where each function is accompanied by a docstring that provides all the necessary information. This information can be easily accessed and displayed in the Python terminal. The following code snippet illustrates this approach:

```

1 import pandas as pd
2 from InstrumentControl.Functions import Coupling_Stability
3 from InstrumentControl.Functions import MZM_transfer
4 from InstrumentControl.Functions import GC_transfer
5 from InstrumentControl.Functions import SG_SA_transfer
6 from InstrumentControl.Functions import SG_SA_Linearing
7 from InstrumentControl.Functions import Laserbeating_RF_Response
8 from InstrumentControl.Functions import Device_characterization
9 from InstrumentControl.Functions import MZM_transfer_Keithley
10

```

```

11 from InstrumentControl.Functions import Voltage_OSA
12 from InstrumentControl.Functions import Sweep_CWF_VNA
13 from InstrumentControl.Functions import Sweep_Power_VNA
14 from InstrumentControl.Functions import Laserbeating_RF_Response_CoBrite
15 from InstrumentControl.Functions import Tobias_MZM_passive_Equalizer
16 from InstrumentControl.Functions import Tobias_MZM_wo_Equalizer
17 from InstrumentControl.Functions import Tobias_MZM_measurement
18 from InstrumentControl.Functions import
    Tobias_MZM_characteristic_in_dependency_of_thermal_biasing
19 import os
20
21
22
23
24 def __init__():
25     print('''
26         ##### Short Function Description #####
27         1) Coupling_Stability:
28             In this function the coupling stability over time will
29             be
30                 measured. The Power Meter is used.
31
32         2) MZM_transfer:
33             In this function the transfer characteristic of a MZM
34             Interferometer
35                 will be measured. For the purpose the MZM will be
36             connected
37                 to a Power Supply that will sweep the Voltage from 0V
38             to a
39                 given voltage. The Output of the MZM will be measured
40             whit
41                 the Power Meter.
42
43         3) MZM_transfer_Keithley:
44             In this function the transfer characteristic of a MZM
45             Interferometer
46                 will be measured. Same as function 2) but with the
47             help of
48                 KEITHLEY Source Meter.
49
50         4) GC_transfer:
51             In this function the LU1000 from Novoptel will be
52             sweepet. The
53                 user will be asked to set a min and max wavelength for
54             the sweep.
55                 The Power Meter will measure the output of the Laser.
56             It can be used
57                 to characterise different optical DUTs.
58
59         5) SG_SA_transfer:
60             In this function the CW Frequency of the Signal
61             Generator will be
62                 sweepet and the Spectrum will be measured with the
63             Spectrum Analyzer.

```



```

87      11) Sweep_Power_VNA:
88          Sweep the Power from the Vectro Analyzer Port by
89          constant continuous
90          waveform(CW) Frequency. Spectrum Analyzer will detect
91          the
92          Signal and save the data in CSV, param and svg files.
93
94      12) Laserbeating_RF_Response_CoBrite:
95          Same as Laserbeating_RF_Response! Will use the CoBrite
96          Laser
97          and the spectrum analyzer. No pewer meter needet.
98
99      13) Tobias_MZM_passive_Equalizer:
100         A Function that is used to measure Tobias Chip(MZM
101         passive Equalizer). The instruments needed are
102         VNA, Power Mether ThorLabs, Source Meter and Power
103         Supply. The Voltage from
104         the Power Supply will be sweeped. At each iteration ,
105         the source meter voltage
106         is swept and for each voltage swept , a voltage
107         measurement is made by the source
108         meter , a power measurement is made by the power meter ,
109         and an S-parameter 4-port
110         measurement is made. The information will be saved in .
111         csv , .svg Files and additional
112             .txt Files with speacial synthax suited for ADS.
113
114      14) Tobias_MZM_wo_Equalizer:
115         A Function that is used to measure Tobias Chip(MZM wo
116         Equalizer). The instruments needed are
117         VNA, Power Mether ThorLabs, Source Meter and Power
118         Supply. The Voltage from
119         the Power Supply will be sweeped. At each iteration ,
120         the source meter voltage

```

```

powermeter.

121             This process is repeated the other way around.
122             Afterwards, the maximum, minimum and the quadrature
123             points are determined and saved,
124             and the MZM is biased to quadrature point.
125             Then, either the voltage or the current of the
126             SystemSourceMeter is swept.
127             Both channels are swept with same voltage or current,
128             and the corresponding current or voltage of both
129             channels is measured.
130             For each voltage or current step, S-parameters are
131             measured by means of the VNA.
132             All measurement results are saved in .snp, .csv, .txt
133             and .svg files.
134             The syntax of the .txt files is adapted to ADS.

135         16) Tobias_MZM_characteristic_in_dependency_of_thermal_biasing:
136             A function which is used to measure an MZM. The
137             instruments needed are
138                 VNA (Anritsu MS4647B with extension), optical
139                 powermeter (ThorLabs S155C),
140                 SystemSourceMeter (Keithley 2612) and 2x Power Supply (
141                 RND Labs 320-KA3005P).
142             The power supplies deliver the voltages fed to the
143             thermal phaseshifters.
144             These voltages are swept and in each step the MZM
145             characteristic is measured via the
146                 the VNA and the SystemSourceMeter.
147             This is repeated for multiple biasing voltages/currents
148             of the regular Si phaseshifters in the MZM.
149             All measurement results are saved in .snp, .csv, .txt
150             and .svg files.
151             The syntax of the .txt files is adapted to ADS.

152         '''
153     )
154
155
156     def MeasFunctions(Instrument ,path):
157         print('',
158             ##### Functions Menu #####
159             0) For short function description = 0
160             1) Coupling_Stability = 1
161             2) MZM_transfer = 2
162             3) MZM_transfer_Keithley = 3
163             4) GC_transfer = 4
164             5) SG_SA_transfer = 5
165             6) SG_SA_Linearing = 6
166             7) Laserbeating_RF_Response = 7
167             8) Device_characterization = 8
168             9) Voltage_OSA = 9
169             10) Sweep_CWF_VNA = 10
170             11) Sweep_Power_VNA = 11
171             12) Laserbeating_RF_Response_CoBrite = 12
172             13) Tobias_MZM_passive_Equalizer = 13

```

```

163     14) Tobias_MZM_wo_Equalizer = 14
164     15) Tobias_MZM_measurement = 15
165     16)
166 Tobias_MZM_characteristic_in_dependency_of_thermal_biasing = 16
167     17) Exit
168
169     ##### Functions Menu #####
170
171     ''')
172
173     FunctNameList = [ '0', '1', '2', '3', '4', '5', '6', '7', '8', '9', '10', '11', '12',
174     , '13', '14', '15', '16', '17' ]
175     FunctName = input( 'Function name: ' )
176     if FunctName in FunctNameList:
177         while FunctName != '17':
178
179             if FunctName == '1':
180                 print('''
181                         Parameters
182
183                         Time : int / float
184                         Time in minutes. How many Min the Power
185                         Detector will measure the
186                         Power.
187                         WaveLength : int
188                         Wavelenght for the Power Meter. It is give in
189                         nm. For example WaveLength = 1310, WaveLength = 1550
190                         path : str
191                         Where to save the file.
192                         name : str
193                         Name of the File. The name is Combinationof the
194                         function name and
195                         name that you give.
196                         For Example: Coupling_Stability_Test.csv
197
198                         Returns
199
200                         None data will be returned. A TXT File with the
201                         Power Detector data will
202                         be created , CSV file whit data will be created and
203                         SVG plot
204                         of the measurment will be created and saved in the
205                         given in 'path' folder.
206                         ''')
207                         print('\n')
208                         Time = float(input('Time = '))
209                         WaveLength = float(input('Power Meter Wavelenth [nm] = '))
210                         name = input('name = ')
211
212                         Coupling_Stability(Instrument ,Time, WaveLength ,path ,name
213
214                     )
215                     break

```

While the function descriptions are not yet complete due to ongoing work on other tasks, they are designed to be explicit enough to provide a clear idea of their intended use. Once you select a function and provide the necessary instrument and measurement parameters in the Python terminal, the information is then passed to Functions.py. In this script, you can find detailed implementations of the measurements. For instance, let's consider an example of a coupling stability measurement. The following code snippet outlines the process:

```

1 import numpy as np
2 import pandas as pd
3 import tkinter as tk
4 from tkinter import filedialog
5 import matplotlib.pyplot as plt
6 plt.rcParams.update({'font.size':22})
7 plt.rcParams["figure.figsize"] = (22,20)
8 import time as t
9 import sys
10 import pyvisa as visa
11 import os
12 from datetime import datetime
13
14
15
16
17 #
=====

18 # Loading Bar
19 #
=====

20
21 def loadingBar(count ,total ,size):
22     percent = float(count)/float(total)*100
23     sys.stdout.write("\r" + str(int(count)).rjust(3,'0')+" /"+
24                             str(int(total)).rjust(3,'0') + '[' + '='*int(percent
25 /10)*
26                             size + ' '*(10-int(percent/10))*size + '] ')
27
28 def loadingBarTwo(count ,total ,size):
29     percent = float(count)/float(total)*100
30     sys.stdout.write("\r" + str(int(count)).rjust(3,'0')+" /"+
31                             str(int(total)).rjust(3,'0') + '[' + '*'*int(percent
32 /10)*
33                             size + ' '*(10-int(percent/10))*size + '] ')
34
35 #
=====

36 # Definitions
37 #
=====
```

```

36
37
38 def Coupling_Stability(Instrument ,Time , WaveLength , path ,name) :
39     """
40
41     Parameters
42     -----
43     Time : int
44         Time in minutes. How many Min the Power Detector will measure the
45         Power.
46     WaveLength : int
47         Wavelenght for the Power Meter. It is give in nm. For example
48         WaveLength = 1310 , WaveLength = 1550
49     path : str
50         Where to save the file .
51     name : str
52         Name of the File . The name is Combinationof the function name and
53         name that you give .
54         For Example: Coupling_Stability_Test.csv
55
56
57     Returns
58     -----
59     None data will be returned . A TXT File with the Power Detector data
60     will
61     be created , CSV file whit data will be created and SVG plot
62     of the measurmend will be created and saved in the given in 'path'
63     folder .
64     """
65
66     for i in range(len(Instrument)):
67         if 'PM100D' in str(Instrument[i]).split('.'):
68             PM = Instrument[i]
69         else:
70             raise Exception(""
71                             The Instrument that you select is not one
72                             of:
73                             — Power Meter
74                             """)
75
76     name = 'Coupling_Stability_' + name
77     nameParam = name + '_Param'
78     nameImage = name + '_Image'
79
80     secs = Time*60
81     TimeVec = np.arange(1,int(secs+1),1)
82     PowerVec = []
83     Headaer = [ 'Time/s' , 'Power/' +PM.ask_PowerUnits() ]
84     for i in range(secs):
85         loadingBar(i,int(secs+1),1)
86         t.sleep(1)
87         PowerVec.append(PM.DefaultPowerMeas(WaveLength))
88
89     PowerVec = np.array(PowerVec , dtype=np.float32)

```

```

86     print('Max Power = ', max(PowerVec))
87     print('Min Power = ', min(PowerVec))
88     data = {Headaer[0]: TimeVec, Headaer[1]: PowerVec}
89     PD = pd.DataFrame(data, columns = Headaer)
90     PD.to_csv(path +'/' + name + '.csv', sep = ',', )
91     Headers, Data, Param = PM.DisplayParamDict('Power')
92     ParamsInst = {}
93     for i in range(len(Param)):
94         ParamsInst[Param[i]] = Data[i]
95     with open(path +'/' + nameParam + '.txt', 'w') as file:
96         for key, value in ParamsInst.items():
97             file.write(key+'\t'+value+'\n')
98
99     plt.figure()
100    plt.plot(PD[Headaer[0]], PD[Headaer[1]])
101    plt.xlabel(Headaer[0])
102    plt.ylabel(Headaer[1])
103    plt.title(name)
104    plt.grid()
105    figManager = plt.get_current_fig_manager()
106    figManager.window.showMaximized()
107    plt.show()
108    plt.savefig(path +'/' + nameImage + ".svg")

```

After completing your measurements, the function will return the measured data in two formats: an .svg image visualizing the data and the raw data saved in either .txt or .csv format and an .txt file with all the instrument settings for the specific measurement. All of this data will be conveniently stored in the directory you specified at the beginning. However, it is imperative to close all the communication ports that were opened between your PC/Laptop and the instruments. Failing to do so may lead to issues. You can achieve this for each instrument in the main.py file using the following code snippet:


```

1 #dont forget to close the instrument connection to the PC. (free the socket
2
3 # To close Anrtisu spectrum Analyzer MS2760A
4 SA.Close()
5
6 # To close Anrtisu Signal Generator MG3694C
7 SG.Close()
8
9 # To close Anrtisu Vectro Analyzer MS4647B
10 VNA.RTL() # Disconnect from VNA so that you can use the VNA afterwards. If
11     not you need to restart the VNA!!!
11 VNA.Close()
12
13 # To close Novoptel Laser LU1000
14 LU.Close()
15
16 # To close Power Meter ThorLabs PM100D
17 PM.Close()
18
19 # To close KEITHLEY Source Meter 2612
20 KA.Close()

```

```
21
22 # To close Power Supply KA3005 All Programmable KA
23 PS.Close()
24 PS1.Close()
25 PS2.Close()
26 PS3.Close()
27
28 # To close Yokogawa Optical Spectrum Analyzer AQ6370D
29 OSA.Close()
30
31 # To close CoBrite Tunable Laser
32 CO.Close()
33
34 # To close AnaPico AG,APPH20G
35 AP.Close()
36
37 # To close 4-Channels Power Suppy GPP4323
38 GPP.Close()
```

In the flow chart 3.1 an example of the Automatic measurements python library construction is again given.



*Abbildung 3.1: Flow Chart of the automatic measurements python scripts*

# 4 Example

In this chapter, we present an example along with a step-by-step picture tutorial on how to utilize the scripts. This example is conducted on the laboratory PC within the laser laboratory, where the Spyder API from the Anaconda distribution is employed.

Let's get started:

- **Step 1: Clone/Download the Automation Measurements Folder from Git**

Begin by cloning or downloading the 'Automation Measurements' folder from the Git repository. This folder contains all the necessary scripts, organized into different folders and sorted by name. In Figure 4.1, you can observe the files available in the 'Automation Measurements' folder, including the essential main.py file.

- **Step 2: Connect the Measurement Instruments and Verify Device Drivers**

Before you begin automated measurements, it's essential to ensure that your measurement instruments are properly connected, and the required device drivers are installed on your system. Follow these steps: 1) Connect your measurement instruments as per your specific setup requirements. 2) Verify that the necessary device drivers for your instruments are correctly installed and recognized by your system. You can check this in the 'Device Manager' on your computer.

(Note: Users should ensure that their instruments are correctly connected and drivers are installed. Specific instructions may vary based on the instruments used.)


In this example, we'll focus on using the Power Meter ThorLabs PM100D as the instrument. Later on, we'll apply the existing Coupling Stability function. The image in Figure 4.2 provides an example of a measurement setup to help you visualize the configuration. With your instruments properly connected and drivers confirmed, you're ready to proceed with the automated measurements

- **Step 3: Open the `main.py` Script**

To begin automated measurements, open the `main.py` script in your Python API environment. In this example, we are using Spyder. You can follow these steps:

1. Launch your Python API (Spyder) on your computer
2. Locate the `main.py` script within the 'Automation Measurements' folder that you previously cloned or downloaded

You can see part of the script's contents in Figure 4.3. Opening 'main.py' is the next step in setting up and running automated measurements. Let's proceed to configure and use the script.



**Abbildung 4.1:** Folders in Auto\_Measurement Git Directory

- **Step 4: Start the Script**

With the `main.py` script open in your Python environment (e.g., Spyder), you're ready to initiate the script sequence by clicking the play button or running the script. Be sure to provide all the required information accurately as prompted. While we've incorporated safety error-catching functions into the scripts, not all scenarios have been thoroughly tested, as I am personally familiar with script operation. If you encounter any errors or can provoke errors that the scripts fail to catch, please report them promptly. Your feedback is invaluable for identifying and resolving any issues to enhance the reliability of the scripts.

- **Step 5: Input the Number of Instruments to Be Used**

Upon clicking the play button to start the script, the first prompt in the Python console



*Abbildung 4.2: Measurement setup*

will ask you to specify how many instruments you intend to use. You will need to enter the correct number accordingly. Figure 4.4 provides a visual representation of what this prompt will look like in the console.

In this example, since i am using the Power Meter ThorLabs PM100D, i will select 1 for one instrument.

- **Step 6: Select the Instrument You Will Use**

After specifying the number of instruments you plan to use, the next prompt in the Python console will ask you to select the instrument from the list of instruments you intend to use. Simply enter the corresponding number for your instrument choice. Figure 4.5 provides a visual representation of what this selection prompt will look like in the console.

After selecting the instrument by entering the corresponding number, the script will proceed to automatically search all the USB and Ethernet ports on your PC for the connected instrument. Once a successful connection is established, you will receive a message indicating that the instrument is connected, such as: **Power Meter ThorLabs PM100D is connected as PM**. Please note the last initials, which, in this case, are PM. With these initials, you can now call functions from the Power Meter ThorLabs PM100D Library directly in the **main.py** file. For instance, if you wish to set the wavelength range on the device, you can easily do so by entering the following command in the Python console: **PM.WaveLength(1550)**. This command will set the wavelength to (1550 nm). This

The screenshot shows the Spyder IDE interface. On the left, the code editor displays `main.py`. The code defines a function `FunctionCall()` which prompts the user for the number of instruments and lists them. The right side shows the IPython console output:

```

Python 3.8.5 (default, Sep 3 2020, 21:29:08) [MSC v.1916 64 bit
(AMD64)]
Type "copyright", "credits" or "license" for more information.
IPython 7.19.0 -- An enhanced Interactive Python.

In [1]: numInst = int(input('Number Instruments = '))
Number Instruments =

```

*Abbildung 4.3: main.py opened in Spyder*

The screenshot shows the Spyder IDE interface. On the left, the code editor displays `main.py`. The code defines a function `FunctionCall()` which prompts the user for the number of instruments and lists them. The right side shows the IPython console output:

```

Python 3.8.5 (default, Sep 3 2020, 21:29:08) [MSC v.1916 64 bit
(AMD64)]
Type "copyright", "credits" or "license" for more information.
IPython 7.19.0 -- An enhanced Interactive Python.

In [1]: runfile('C:/Users/SCT-Labor/Desktop/Python Automation skips/main.py', wdir='C:/Users/SCT-Labor/Desktop/Python Automation skips')
Number Instruments =

```

*Abbildung 4.4: How many instruments will be used*

The screenshot shows the Spyder Python IDE interface with the following details:

- File Menu:** File, Edit, Search, Source, Run, Debug, Console, Projects, Tools, View, Help.
- Toolbar:** Includes icons for Open, Save, Run, Stop, and Help.
- Code Editor:** The main area displays Python code in `main.py`. It includes comments like `##### Select which instruments will be used #####` and `Table of instruments` followed by a list of connected equipment.
- Console:** Shows the output of the code execution, including the table of instruments and the result of the `numInstr = int(input('Number Instruments = '))` command.
- Output Window:** Displays the table of instruments and the selected instrument name.
- Help:** A context help box is open over the code, explaining how to get help for objects and how to enable automatic help for objects.
- Bottom Status Bar:** Shows "LSP Python: ready", "conda: base (Python 3.8.5)", "Line 59, Col 6", "UTF-8", "CR/LF", "RW", and "Mem 56%".

**Abbildung 4.5:** Which instruments will be used. Select from list

This step simplifies the process of connecting to instruments and accessing their functions, allowing you to control and configure instrument settings with ease.


- Step 7: Choose the Location to Save Measurement Data

After connecting to the instruments and configuring your measurements, the next important step is choosing where on your PC you want to save the extracted measurement information. As mentioned in Chapter 1, the script will generate three types of files for each measurement:

1. An .svg picture visualizing the measurement.
  2. An .txt file containing all the parameters you set for the measurement.
  3. A .csv file containing all the measured data, allowing you to plot the measurement graph as needed.

To facilitate this process, an automatic pop-up window (as shown in Figure 4.6) will appear, prompting you to navigate to the folder where you want to save this valuable data. This step ensures that your measurement data is efficiently organized and saved to the location of your choice for future analysis and reference.

You can now write your own measurement function or you can use one of the pre-programmed functions. To do so uncomment the following line in the `main.py` file.



**Abbildung 4.6:** Select folder on your PC to save the data

```

1 #Main Function. Instrument is a list of selected instruments for the
   measurement
2 MeasFunctions(Instrument ,path)

```

- **Step 8: Use Pre-Programmed Measurement Functions (Optional)**

If you prefer to use the already pre-programmed measurement functions, you can do so with ease. Simply use the `MeasFunctions(Instrument, path)` function in the script. However, if you wish to skip this step and customize the script further, you have the flexibility to comment out this line based on your preference. For example, in my example scenario, I will use the `MeasFunctions(Instrument, path)` function to perform the existing 'coupling\_stability' measurement. Upon using the `MeasFunctions(Instrument, path)`, a new set of questions will appear in the Python console. The first question will prompt you to choose from one of the already existing functions, as shown in Figure 4.7.

- **Step 9: Set Measurement Data Parameters**

After selecting the pre-programmed measurement function, the script will prompt you to set all the necessary data parameters for the measurement. The required parameters for your chosen function will be clearly outlined in the docstring, which will be printed in the Python console. This docstring provides detailed guidance on what data needs to be configured for the measurement. For visual reference, Figure 4.8 illustrates how this

The screenshot shows the Spyder Python IDE interface. The top menu bar includes File, Edit, Search, Source, Run, Debug, Consoles, Projects, Tools, View, Help. The toolbar has icons for file operations like Open, Save, Run, Stop, and Help. The left sidebar shows project files: main.py, FunctionCall.py, and Functions.py. The main code editor window contains Python code for instrument selection and connection. The right side features a 'Source' tab, a 'Console' tab with output, and a 'Usage' help panel. The console output shows the execution of the script, including instrument names and their connections.

```
22
23 numInstr = int(input('Number Instruments = '))
24 print("***** How many instruments will be used *****")
25
26
27
28
29
30
31
32
33
34
35 print("***** Select which instruments will be used *****")
36
37 You will be asked to give the name of the instruments
38 manually. For example Instrument name: MS4647B for
39 Anritsu Signal Analyzer.
40
41 Table of Instruments
42 Anritsu spectrum Analyzer MS2760A = 1
43 Anritsu Signal Generator MG3694C = 2
44 Anritsu Vetro Analyzer MS4647B = 3
45 Power Meter ThorLabs PM100D = 4
46 Novoptel Laser LU1880 = 5
47 Yokogawa Optical Spectrum Analyzer AQ6378D = 6
48 KEITHLEY Source Meter 2612 = 7
49 Power Meter PM100 = 8
50 CoBrite Tunable Laser = 9
51 Anpico AG_APPHZB6 = 10
52 4-Channels Power Supply GPP4323 = 11
53
54
55
56
57
58
59
60
61
62
63
64
65
66
67
68
69
70
71
72
73
74
75
76
77
78
79
80
81
82
83
84
85
86
```

```
InstrumentNameList = ['1','2','3','4','5','6','7','8','9','10','11']
Instrument = []
for i in range(numInstr):
    instrumentName = input('Instrument name: ')
    if instrumentName in InstrumentNameList:
        Instrument.append(instrumentName)
    else:
        raise ValueError('Invalid Instrument Selected')

for elements in listInstruments:
    if elements == '1':
        SVA = InstInit(elements)
        Instrument.append(SVA)
        print('Anritsu Spectrum Analyzer MS2760A is connected as SA')
    elif elements == '2':
        SG = InstInit(elements)
        Instrument.append(SG)
        print('Anritsu Signal Generator MG3694C is connected as SG')
    elif elements == '3':
        VNA = InstInit(elements)
        Instrument.append(VNA)
        print('Anritsu Vetro Network Analyzer MS4647B is connected as VNA')
    elif elements == '4':
        PM = InstInit(elements)
```

Usage

Here you can get help of any object by pressing **Ctrl+I** in front of it, either on the Editor or the Console.

Help can also be shown automatically after writing a left parenthesis next to an object. You can activate this behavior in **Preferences > Help**.

Variable explorer Help Plots Files

Console 1/A Console 4/A

```
4-Channel Power Supply GPP4323 = 11

#####
Select which instruments will be used #####
Instrument name: 4
Thorlabs,PM100D,P0024970,2,6,0
Thorlabs,PM100D,P0024970,2,6,0
Power Meter ThorLabs PM100D is connected as PM

#####
Functions Menu #####
0) For short function description = 0
1) M2M_characteristic = 1
2) M2M_transfer = 2
3) M2M_Transfer_Keithley = 3
4) M2M_Transfer = 4
5) SG_SpectrumAnalyzer = 5
6) SG_SA_linearize = 6
7) Laserbeating_RF_Response = 7
9) Tobias_M2M_characterization = 8
9) Voltage_OSA = 9
10) Sweep_VNA = 10
11) Sweep_Power_Meter = 11
12) Tobias_M2M_RF_Response_CoBrite = 12
13) Tobias_M2M_Equalizer = 13
14) Tobias_M2M_Wo_Equalizer = 14
15) Tobias_M2M_measurement = 15
16)
Tobias_M2M_characteristic_in_dependency_of_thermal_biasing = 16
17) Tobias_M2M_V_T_pi_measurement = 17
18) Tobias_M2M_measurement_4_channel_voltage_source =
19) Tobias_thermal_phaseshifter_characterization = 19
20) Tobias_plasmonic_phaseshifter_DC_characterization =
21) Exit

#####
Functions Menu #####
Function name: |
```

**Abbildung 4.7:** Select one of the pre-programmed measuring functions

process appears in the Python console. This step ensures that you have full control over configuring the measurement according to your specific requirements, with the docstring serving as a comprehensive guide.

- Step 10: Commence the Measurement

Once you have provided all the necessary information and configured the measurement parameters, the measurement process will commence. In this example, it involves measuring the power at a given wavelength for a set duration. You will be able to observe the progress of the measurement in real-time. Figure 4.9 provides a visual representation of what the measurement process looks like in the Python console. You will notice a loading bar indicating the progress of the measurement, along with a countdown showing the number of sweeps completed and the remaining sweeps to be performed. This step marks the execution of the measurement, allowing you to track its progress and gather valuable data for analysis.

- Step 11: Close Open Sockets and Disconnect Instruments

Upon completing your measurements, it is crucial to close all open sockets. These sockets act as communication channels between your PC and the instruments. Failure to close these sockets can lead to issues, particularly with instruments like the VNA (Vector Network Analyzer), which may require a restart if the socket is not released. For the VNA,

The screenshot shows a Python development environment with several windows open. On the left, a code editor displays a script named `FunctionCall.py` containing Python code for instrument selection and measurement setup. The code includes a table of instruments and logic to select which ones will be used based on user input. On the right, a help window titled "Usage" provides instructions on how to use the `help` command. Below it, a console window shows the output of running the script, including the table of instruments and the selection process. The console also shows the user interacting with the script to set measurement parameters like time, wavelength, and file name.

```

File Edit Search Source Run Debug Consoles Projects Tools View Help
C:\Users\SCTLabor\Desktop\Python Automation scripts\Functions.py
Source Console Object
Usage
Here you can get help of any object by pressing Ctrl+I in front of it, either on the Editor or the Console.
Help can also be shown automatically after writing a left parenthesis next to an object. You can activate this behavior in Preferences > Help.
Variable explorer Help Paths Files
Console 1/A Console 4/A
14 Tobias_M2M_characterization = 16
15 Tobias_M2M_in_dependency_of_thermal_biasing = 16
16 Tobias_M2M_V_pi_measurement = 17
17 Tobias_M2M_measurement_4_channel_voltage_source =
18 18 Tobias_thermal_phaseshifter_characterization = 19
19 Tobias_plasmonic_phaseshifter_DC_characterization =
20 20 Tobias_M2M_measurement_4_channel_voltage_source =
21 Exit
#####
Functions Menu #####
Function name: 1
Parameters
Time : int,float
Detector will measure the
Power
Wavelength : int
Wavelength for the Power Meter. It is
give in nm. For example Wavelength = 1310, Wavelength = 1550
path : str
where to save the file.
name : str
Name of the File. The name is
Combinationof the function name and the
name that you give.
For Example: Coupling_Stability_Test.csv
Returns
None
None data will be returned. A TXT File with
the Power Detector data will
be created, CSV file whit data will be
created and SVG plot
of the measurment will be created and saved
in the given in 'path' folder.
Time = |
Python console History
LSP Python ready conda: base (Python 3.8.5) Line 99, Col 6 UTF-8 CR/LF R/W Mem 51%

```

**Abbildung 4.8:** Set the required measurement parameters. In this section I am asked to set time, wavelength and name for my measurement files

this means recalibrating the instrument, which can be time-consuming. To close the socket associated with a specific instrument, you can use a command in the following format: `InstrumentName.Close()`. Replace 'InstrumentName' with the initials provided by the Python console when connecting the instrument, followed by 'Close'. For example, to close the socket for the Power Meter ThorLabs PM100D, you can enter `PM.Close()`. Figure 4.10 provides a visual reference on how to close a socket and disconnect an instrument in the Python console. This step ensures that all communication channels are properly closed, preventing potential issues and ensuring a smooth workflow for future measurements.

Once your measurement is complete, you can conveniently access your measurement results in the folder you specified at the beginning of the process (see Figure 4.11 for reference). The data generated and saved includes:

1. An .svg file that visualizes the measurement
2. An .txt file containing all the parameters you configured for the measurement
3. A .csv file that holds all the measured data, enabling you to plot and analyze the results as needed (see Figure 4.12)

This step enables you to access and work with your measurement data for further analysis and interpretation.

```

 115         elif data == SerialNum[1]:
 116             PS2 = PS
 117             print('Power Supply_KA3005P is connected as PS2')
 118             Instrument.append(PS2)
 119             if data == SerialNum[2]:
 120                 PS3 = PS
 121                 print('Power Supply_KA3005 is connected as PS3')
 122                 Instrument.append(PS3)
 123
 124
 125     elif elements == 'C':
 126         C0 = InstInstruments()
 127         Instrument.append(C0)
 128         print('CoDrite Tunable Laser is connected as C0')
 129     elif elements == 'AP':
 130         AP = InstInstruments()
 131         Instrument.append(AP)
 132         print('Anapico_APPH200 is connected as AP')
 133     elif elements == 'GPP':
 134         GPP = InstInstruments()
 135         Instrument.append(GPP)
 136         print('GW InsterK_GPP4523 is connected as GPP')
 137     else:
 138         raise ValueError('Invalid Instrument Selected')
 139
 140
 141 # =====#
 142 # Select Path for Save Data (pop up window)
 143 # =====#
 144 root = tk.Tk()
 145 path = filedialog.askdirectory(parent = root)
 146 root.destroy()
 147
 148
 149
 150
 151 # =====#
 152 # Functions
 153 # =====#
 154 # Main Function:Instrument is list of selected instruments for the measurement
 155 MeasFunctions(Instrument,path)
 156
 157
 158
 159 #don't forget to close the instrument connection to the PC. (free the socket)
 160
 161 # ## # SA.Close()
 162 # ## # SG.Close()
 163 # ## # LU.Close()
 164
 165 # ## # PM.Close()
 166 # ## # KA.Close()
 167 # ## # GPP.Close()
 168
 169 # ## # VNA.RTL()
 170 # ## # VNA.Close()
 171
 172 # ## # PS.Close()
 173 # ## # PS1.Close()
 174 # ## # PS2.Close()
 175 # ## # PS3.Close()
 176 # ## # CO.Close()
 177 # ## # AP.Close()
 178 # ## # GPP.Close()
 179 # ## # GPP.Close()

```

The screenshot shows the Spyder Python IDE interface. On the left, the code editor displays a script named 'main.py' with several sections of code. A help pop-up window titled 'Usage' is open in the top right corner, providing information on how to get help for objects in the editor or console. On the right, the 'Console' tab is active, showing command-line input and output. The output includes lines like 'Tobias\_thermal\_phaseshifter\_characterization = 19', 'Tobias\_plasmonic\_phaseshifter\_DC\_characterization = 20', and 'Exit'. Below the console, a 'Function name: 1' section shows detailed documentation for a function, including parameters, detector measurements, and returns. The bottom right of the interface shows tabs for 'Python console' and 'History'.

**Abbildung 4.9:** The measurement is set and running. The user can see how much is left and how many are done through the load bar function

In this example, we've demonstrated the streamlined process of using the Automated Measurements Python script to perform measurements with ease. While we've used only one instrument for this illustration, it's important to note that the script supports multiple instruments, allowing you to harness the full range of capabilities if needed. Moreover, we've showcased the convenience of pre-programmed measurement functions. These functions are ready to use for quick measurements. However, if you have specific measurement needs or wish to create custom functions, you have the flexibility to reach out for assistance, or if you're proficient in Python, you can build custom functions yourself. The Automated Measurements Python script empowers you to efficiently control and analyze measurements, making it a valuable tool in various laboratory and research settings. Feel free to explore, adapt, and customize the script to suit your unique requirements. Your measurement capabilities are now in your hands, whether you need swift measurements or tailor-made solutions.

Should you have any questions or require assistance with custom functions, don't hesitate to reach out. Happy measuring!

```

115 elif data == SerialNum[1]:
116     PS2 = PS
117     print('Power Supply, K43085P is connected as PS2')
118     Instrument.append(PS2)
119     if data == SerialNum[2]:
120         PS3 = PS
121         print('Power Supply K43085 is connected as PS3')
122         Instrument.append(PS3)
123
124 elif elements == 'C0':
125     C0 = InstInstruments
126     Instrument.append(C0)
127     print('CoDrite Tunable Laser is connected as C0')
128 elif elements == 'AP':
129     AP = InstInstruments
130     Instrument.append(AP)
131     print('AmaPico AG_APH280 is connected as AP')
132 elif elements == 'GRP':
133     GRP = InstInstruments
134     Instrument.append(GRP)
135     print('GW InsteK_GPP4323 is connected as GPP')
136 else:
137     raise ValueError('Invalid Instrument Selected')
138
139 # =====
140 # Select Path for Save Data (pop up window)
141 # =====
142 root = tk.Tk()
143 path = filedialog.askdirectory(parent = root)
144 root.destroy()
145
146
147
148
149
150
151 # =====
152 # Functions
153 # =====
154
155 def MainFunctionInstrument(instrument, path):
156
157
158     #don't forget to close the instrument connection to the PC. (free the socket)
159
160     # # # # SA.Close()
161     # # # # SG.Close()
162     # # # # LU.Close()
163     # # # # Close()
164
165     PW.Close()
166     # KA.Close()
167     # GPP.Close()
168
169     # VNA.RTL()
170     # VNA.Close()
171
172     # PS.Close()
173     # # PS1.Close()
174     # PS2.Close()
175
176     # # # # Close()
177     # # # # CO.Close()
178     # # # # AP.Close()
179     # # # # GPP.Close()

```


Abbildung 4.10: Disconnect the device from the mains after use

```

115 elif data == SerialNum[1]:
116     # # # # Measurements_1_10_23
117     # Datei Start Freigeben Ansicht
118     # Schnellzugriff Kopieren Einfügen Anhängen Zwischenablage Verschieben Kopieren nach... Verknüpfung einfügen Löschen Umbenennen Neuer Ordner Organisieren Neu Einfacher Zugriff Eigenschaften Bearbeiten Verlauf Alle auswählen Nichts auswählen Auswahl umkehren Auswählen
119     # # # # Desktop
120     # # # # Downloads
121     # # # # Dokumente
122     # # # # Bilder
123     # # # # Dieser PC
124     # # # # Pictures_Docu
125     # # # # SVG_V_th_sweep
126     # # # # SVG_V_th_sweep
127     # # # # Tabular_idc_voltages_currents_powers_
128
129 elif
130
131 else:
132     # =====
133     # Select
134     # =====
135     root = t
136     path = des
137
138
139
140
141
142     # Select
143     # =====
144     root = t
145     path = des
146
147
148
149
150
151     # =====
152     # Functions
153     # =====
154
155
156     # Main Function
157     # =====
158     # don't forget to close the instrument connection to the PC. (free the socket)
159
160     # # # # SA.Close()
161     # # # # SG.Close()
162     # # # # LU.Close()
163     # # # # Close()
164
165     PW.Close()
166     # KA.Close()
167     # GPP.Close()
168
169     # VNA.RTL()
170     # VNA.Close()
171
172     # PS.Close()
173     # # PS1.Close()
174     # PS2.Close()
175
176     # # # # Close()
177     # # # # CO.Close()
178     # # # # AP.Close()
179     # # # # GPP.Close()

```

Abbildung 4.11: Data saved on PC/Laptop after measurement



**Abbildung 4.12:** Contents of the .csv file